

Automation-Special machine construction

Krommebeekstraat 24 – BE 8930 MENEN Tel +32 56 52 14 80 - Fax +32 56 53 00 15 www.vaskon.com

## Deepwell shaker

(shaking machine for DNA research)





Pict. 1: cam box with package of products on top of it that can be shaken. The complete unit is protected by a locked door made of polycarbonate.



Pict. 2: a package of products is also put at the bottom.



Pict. 3: detail of products, shaken perfectly linear (test tubes containing agricultural seeds)  $(amplitude\ 4cm-frequencu\ adjustable\ from\ 0\ to\ 30\ Hz).$ 



Pict. 4: Elektrical desk with push buttons:

- main switch with indication light
- emergency stop
- rotary button for frequency ajdustment
- start- and stop button
  adjustable timer with LCD-display (upper side)

Drive : asynchrone motor (speed frequency adjustable through frequency regulator and manual

rotary button)

: through a mechanical cam Transmission

Amplitude : 40 mm

Frequency : adjustable 0-30 Hz

: adjustable through LCD-display Time duration

Connection : 230 V single-phase

Krommebeekstraat 24 BE - 8930 MENEN Tel +32 56 52 14 80 Fax +32 56 53 00 15 www.vaskon.com